

"WATERDROP DIAMOND FILL"

Luxury care that fills the skin with excellent moisturizing power and diamond luster.

DIASONO·310



Waterdrop Diamond Fill

DIASONO•310

SONO•310 's 3/10MHz **Modulated Ultrasonic Wave** and **Diacell_{peel}** 's **Micro DermaBration** feature create a perfect combination to provide safe and effective treatment and quick recuperative power.



Waterdrop Diamond Fill

DIASONO•310

PROGRAM NAME

“Waterdrop Diamond Fill”

Luxury care that fills the skin with excellent moisturizing power and diamond luster.



FUNCTIONS

3/10Mhz multi ultrasound
DIACELL Peel
Electroporation(E/P)



PROGRAM CONCEPT

- Add features to further enhance the efficacy of 3/10Mhz multi ultrasound and improve the dullness of skin care practices.
- Increase the efficacy and customer satisfaction by simply changing the tips of physical peeling and E/P.




RESULTS

Recovery, regeneration, lifting, acne, tone-up, calming, swelling & bruising relief, pigmentation improvement and wrinkle improvement




SONO·310

3MHz Ultrasound




With 3 million vibrations per second, transmits micro-massage effect up to 10mm depth of skin

10MHz Ultrasound



With 10 million vibrations per second, transmits micro-massage effect up to 3mm depth of skin

3/10MHz Ultrasound



Multi frequency ultrasound, activate cells and extracellular materials

SONO·310 Special function

Cellular and extracellular substrate activation

Utilize the pressure caused by the difference between two different frequencies of ultrasound to provide a massage effect for the skin

01

1,000th of a second (millisecond) cross method

One probe can alternately use two different frequencies of ultrasound, so that the heat can reach the deeper layers of the skin

02

Promote blood circulation

Intercellular motor friction expands blood vessels, increasing blood circulation.

03

Increased oxygen supply

Create intercellular microflow to increase oxygen supply

04

promote the production of collagen

10MHz nano-cell massage can prevent the occurrence of excessive cross-link that causes the skin to lose elasticity, and at the same time promote the production of collagen

05

MMPs ↓

[Proteolysis enzyme]

- ☑ Inflammatory inhibitor
- ☑ Treatment for chronic diseases
- ☑ Scar prevention

MMPs Active Control (decrease)

Proteolysis enzyme

Enzymes that break down functioning proteins+ regulate the activation of new proteins and cells.

It is effective for sensitive and inflammatory skin.

GAGs ↑

[Pethydroxy polysaccharide]

- ☑ Elasticity enhancement
- ☑ Material metabolism enhancement
- ☑ Skin Moisturization

GAGs Content Increase

Moisturizing factor

Combining with moisture in the skin tissue, enhances the metabolism and moisturizes the skin

Helps to have healthy and elastic skin.

HSPs ↑

[thermal shock protein]

- ☑ Skin immunity improvement
- ☑ Regeneration Capacity Improvement
- ☑ pigmentation, wrinkle prevention

HSPs amplification and activation

Thermal shock protein

In the event of skin failure, it is secreted for quick recovery, improving the skin's immune and regenerative ability.

Increases up to 1800% to help repair damaged skin.

Diacell peel

Physical peeling created with tornado tip + medical scaling with AHA and BHA. 3/10Mhz multi ultrasound treatment that brings synergy effects

01
Tornado diamond tip

- Physical peeling
- Exfoliation
- Cell Regeneration

02
Electroporation Tip


- Drug penetration

03
Aqua Peel Tip

- excess sebum removal
- removal of follicles
- Remove makeup residue




DIASONO•310 has a HIGHER ultrasonic penetration!




Before treatment

Ultrasonic penetration rate is low.



Step 1

Exfoliation with diacell peel



Step 2

Increases the penetration rate by applying 3/10Mhz ultrasound. Activate the extracellular materials.

Diacell peel The effect on the skin



Ease of drug penetration

Effectively remove dead keratinocytes, promote regeneration and penetration of the solution



Lifting

Rearrange cells to make the dermis healthy, effectively alleviate fine lines and reduce pores



Whitening

When the keratins fall off, the melanin will fall off together with a whitening effect

DIASONO·310 luxury package care



STEP.1

Diacell^{peel}

- It is used as the first step of various procedures to prevent the keratin from causing the laser to not fully exert its effect
- Utilizing the tornado diamond head and the suction technology of wet compress to safely remove dead keratins, acting on the dermis, and promoting collagen production



STEP.2

E/P

- After the Diacell peel operation is completed, the E/P function can allow the corresponding ampoule to penetrate into the deep layer of the skin in large quantities to improve various skin problems
- Apply electrical energy to the skin to temporarily open the skin, allowing the effective ampoule to penetrate.



STEP.3


SONO·310

- There are two different high-frequency ultrasounds, which can act on different skin levels, promote cell activity, stimulate a variety of biochemical reactions, and return the skin to the state of baby skin
- Coupled with the combined effect of the first and second steps, you can feel the amazing improvement of the skin



Before & After

DIASONO-310 Clinical Photo





THROUGH CONTINUOUS
RESEARCH AND DEVELOPMENT,
MEDICAL BEAUTY DEVICES ARE
PROVIDED TO KEEP UP
WITH THE LATEST TRENDS.

SPECIFICATION

Product Name	DIASONO · 310 ULTRASOUND
Type / Model	LY-DIASONO-310
Size / Weight	275(Weidth)*291(Height)*350(Length) / 9Kg
Power Supply	100~240Vac, 50/60Hz, Power Consumption 80Wmax
Suction Pressure	MAX 600 mmHg Suction LEVEL : 10LEVEL Flow LEVEL : 10LEVEL
E/P	LEVEL : 10LEVEL LED : BLUE LED / LED ON When Running
Ultrasound	Frequency Selection 3MHz, 10MHz, 3/10MHz LEVEL : 10LEVEL Running Tiome : MAX 60min MODE : 8 Program Mode / 4 User Mode LED : BLUE LED / LED ON when Running

MEDICAL SKIN CARE
Lyzerr

Lyzerr Co., Ltd #401, 4F, dongseoul bldg, 466-21, gunja-dong, gwangjin-gu, seoul, korea
TEL : +82-70-4801-6860 FAX: +82-2-2038-4980 E-MAIL: info@lyzer.co.kr

*puREyoung is Lyzerr's beauty device brand.